Ecusit-Composite Repair

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English

Product description

Ecusit-Composite Repair is a light-cure bonding agent especially developed for the repair and shade correction of Composite restorations. The material allows a minimally invasive and time-saving procedure for direct intraoral use. Since Ecusit-Composite Repair was formulated without methyl methacrylate, it is odor-free and exhibits excellent handling properties.

Indications

Shade correction of discolored, caries-free Composite restorations, repair of damaged or abraded Composites.

Contra-indications

Do not use directly on exposed pulp or in cases of allergies to any of the ingredients.

Basic Safety Instructions

- ▶ For dental use only!
- ▶ Keep out of reach of children!
- ▶ Avoid contact with the skin! In the event of inadvertent skin contact, wash the affected area thoroughly with soap and water.
- ▶ Avoid contact with eyes! In the event of inadvertent contact with eyes, immediately rinse the eyes with plenty of water and consult a physician if necessary.
- ▶ Do not apply material directly onto exposed dentine surfaces.

Side effects

There are no known side effects to date. Isolated cases of hypersensitivity to components of the material cannot be excluded.

Instructions for use

- ▶ Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be a minimum of 400mW/cm². Place the light as close as possible to the material.
- ▶ Observe the manufacturer's instructions for any other materials that are used with Ecusit-Composite Repair.

Recommended Use

1. Roughen surface of the Composite to be repaired with a coarse grit disc or a diamond burr (peak-to-valley roughness height: approx. 150 µm). Smooth fractured margins.
2. Position the rubber dam.
3. Etch exposed enamel.
4. Apply a suitable bonding agent to exposed dentine. Observe the manufacturer's instructions for the bonding agent used. If necessary, etch the dentine.
5. Apply one drop Ecusit-Composite Repair with a brush to the area in need of repair and work in thoroughly for at least 30 seconds. In the case of restorations of chemically cured Composites, it is recommended to work in the material for approximately 60 seconds.

Note: Strict adherence to the time indicated is necessary since this is the time the material needs for sufficient penetration of the Composite.

6. Blow the Ecusit-Composite Repair with a gentle air current and cure for at least 20 seconds.
7. Apply, cure and finish the Composite restorative material according to manufacturer's instructions.

Note: The first Composite layer must be applied directly to the dispersion layer using firm pressure. DMG recommends the use of light-curing Composites such as EcuSphere.

Storage and shelf life

- ▶ Store at room temperature (15–25 °C/59–77 °F) in a dry place.
- ▶ Close container immediately after use to avoid curing of the material by ambient light.
- ▶ Do not use after the expiration date.

Composition

Bis-GMA-based dental resins, additives, catalyst. Contains no methyl methacrylate.

Presentation form

1 bottle @ 5 ml bond, 10 brushes, 1 mixing palette	REF 210873
--	------------

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